Installing Python

If Python is already installed on your system, you can skip to the next section.

Step 1: Install Homebrew

Homebrew is the easiest way to install Python on macOS. If you already have Homebrew installed, skip this step.

Press Command+Space on your keyboard. This should open Spotlight and then type in Terminal and press enter.

Paste the following command into the terminal and press Enter:

/bin/bash -c "\$(curl -fsSLhttps://raw.githubusercontent.com/Homebrew/install/HEAD/install.sh)

Step 2: Install Python via Homebrew

In the same terminal, paste the following command and press Enter:

brew install python

Run the command below to verify Python is installed correctly (This should display a version like Python 3.##.#) If nothing appears, try restarting your computer and running it again.

python3 --version

Step 3: Upgrade pip

Paste the command below in the terminal and press Enter:

python3 -m pip install --upgrade pip

Setting up Virtual Environment

- 1. Download the zip file from Google Drive.
- Extract the contents and save them to your Desktop.
- 3. In the terminal, paste the command below and press enter cd ~/Desktop/aln-tool/aln-tool
- 4. Open the terminal and paste the following command below line by line and press enter.

python3 -m venv seq-align

source ./seq-align/bin/activate

cd ./BinghamCanyon

pip install -r requirements.txt


python3 display.py

A desktop browser should appear.

Running the Tool

Q TRY TO AUTOMATICALLY FIND EXE_PATH

- 1. Press the Button:
 - You should see a green alert box saw that your path was set successfully.
 - i. If not and instead see a "Failed to Detect" alert, please press "Choose File" button which corresponds to 1 in the image below. You will be prompted with a file explorer, and select the folders below in this order (Note that specifically ChimeraX might be named differently in the applications folder such as ChimeraX-1.9-rc2024.12.03, so please choose the one you have): Applications ChimeraX- Contents MacOS ChimeraX. Do the same for button 2 navigating: Applications Jalview.app Contents MacOS JavaApplicationStub



2. Press the Button:



- 3. To run the button Press
 - a. You should see a loading icon. This might take ~ 8 minutes, so please wait. You will know the tool is complete when jalview and ChimeraX opens and the loading icon is gone.

Things to Note

- 1. After Changing any of the settings in Method Settings, please press Save Settings, or else the program won't use them.
- 2. Currently, only 1 job can be ran at a time
- 3. To restart the program, you can simple refresh the browser. To do a hard restart go to the terminal and press ctrl-c twice. Then run python3 display.py